



## Return Cartridge Instructions


Eclipse™ CO<sub>2</sub> Adsorbent cartridges are recyclable and are to be returned to the manufacturer for reprocessing. Exhausted cartridges, those having been exposed to CO<sub>2</sub> during use, are now lithium carbonate and do not need to be shipped as corrosive materials.

If using the original shipping box place, fold and place this return sheet (with return address visible) into clear packing list envelop (supplied) and place over the box label shown below.



**FOLD LINE**

**DO NOT DISPOSE OF CARTRIDGE - THE ADSORBENT MATERIAL IS RECYCLABLE**



**USE TO RETURN EXHAUSTED CARTRIDGES TO**

**MICROPORE INC.**

**1000 Konica Drive  
Elkton, MD 21921**

[www.EclipseCO2Adsorbent.com](http://www.EclipseCO2Adsorbent.com)